




Medication Misuse Prevention Solution for the Visually Impaired

App Prerequisites




Android: TalkBack

- How to Enable TalkBack : Settings > Accessibility > Enable TalkBack



iOS: VoiceOver


- How to Enable VoiceOver : Settings > Accessibility > Enable VoiceOver




Medication Misuse Prevention Solution for the Visually Impaired

How to use the Allimi App Features


Allimi App 1 - Home page



Allimi App 2 - Medicine Info page




Allimi App 3 – Medication Detail Info page




Detail Info


Check various drug information such as drug images, efficacy, and dosage.



Allimi App 4-1 - Records Page(Select date)






Allimi App 4-2 - Records Page(Intake Status)



Records Info

Change whether or not you are taking the medication and delete your medication history.




Medication Summary Info TTS

Drug name
ADVIL COLD AND SINUS

The Food and medication Administration classification
Pain reliever/fever reducer - Nasal decongest...

Allimi App 5 – Settings page



Bedtime Settings
Set a blackout time to prevent medication from being scheduled (so the alarm doesn't go off at bedtime)


Settings Language
Multilingual available (Korean, English)



Medication Misuse Prevention Solution for the Visually Impaired

How to use the Allimi App

How to Allimi App 1



RFID Tag Scan

Scan the RFID tag attached to the medicine with the NFC antenna on the back of your phone.


How to Allimi App 2



Bring drug information

The drug information notification app recognizes RFID tags and brings drug information.

How to Allimi App 3




This medicine is ADVIL COLD AND SINUS.
The class is Pain reliever/fever reducer
- Nasal decongestant.

Voice guidance for drug summary information

Summarized drug information provides audible information about the drug name and the Drug Safety classification stored on the RFID tag. The summarized Drug Information page allows you to select abbreviated drug information, detailed information, and whether to take the medication.

How to Allimi App 4




Voice guidance for detailed information on medicines


You can confirm various drug information, including drug images, processing effects, dosage, and precautions.

How to Allimi App 5

[Medication Records]




[Completed Dosage]




- If taken immediately
- If taken after the scheduled dose

[Scheduled to take]



- If you have made a reservation to take

[Missing a Dosage]



- If you miss a dose within 10 minutes of your scheduled appointment, it will be marked as missed.

Medication History

When you click on the dosage record, you can hear the information about canceled dosage (when dosage is completed), completed dosage (when scheduled to be taken, when dosage is missed), and summary information again.